

# Archives systems

## for blocks and slides

### Technical features

**Mobile storage cabinet for blocks** HL-5402-14

Number of trays: 336, 28 double drawers/column

Rows 12, structure with wheels

Capacity (for singol drawer) about 55 blocks


Capacity (for unit) about 36,960 blocks

Unit dimensions (WxDxH) 107 x 56 x 194cm

Single unit weight 190 Kg

Single unit weight (fully loaded) about 480 Kg

Structure in steel white and colour drawers on request



**Mobile storage cabinet for slides** HL-5401-13

Number of drawers 312 trays

Rows 12, structure with wheels

Capacity (for drawer) about 400 slides


Capacity (for unit) about 125,000 slides

Unit dimensions (WxDxH) 107 x 56 x 194cm

Single unit weight 190 Kg

Single unit weight (fully loaded) about 680 Kg

Structure in steel white and colour drawers on request



**Static storage cabinet for blocks** HL-6602-14

Number of trays: 336, 28 double drawers/column

Rows 12

Capacity (for singol drawer) about 55 blocks

Capacity (for unit) about 36,960 blocks


Unit dimensions (WxDxH) 107 x 56 x 184cm

Single unit weight 190 Kg

Single unit weight (fully loaded) about 480 Kg

Structure in steel white and colour drawers on request

Transpallet use



**Static storage cabinet for slides** HL-6501-13

Number of trays: 312 trays

Rows 12

Capacity (for drawer) about 400 slides

Capacity (for unit) about 125,000 slides


Unit dimensions (WxDxH) 107 x 56 x 184cm

Single unit weight 190 Kg

Single unit weight (fully loaded) about 680 Kg

Structure in steel white and colour drawers on request

Transpallet use



**Storage cabinet for blocks - Half capacity** HL-6602-04

Number of trays: 112, 28 double drawers/column

Rows 4

Capacity (for singol drawer) about 55 blocks


Capacity (for unit) about 12,320 blocks

Unit dimensions (WxDxH) 107 x 55 x 67cm

Single unit weight 80 Kg

Single unit weight (fully loaded) about 180 Kg

Colour structure White, colour drawers on request



**Storage cabinet for slides - Half capacity** HL-6501-04

Number of trays: 104

Rows 4

Capacity (for singol drawer) about 400 slides


Capacity (for unit) about 41,600 slides

Unit dimensions (WxDxH) 107 x 55 x 67cm

Single unit weight 80 Kg

Single unit weight (fully loaded) about 250 Kg

Colour structure White, colour drawers on request




Optional for all version: locking doors at request.

Distributed by:




## Archives systems for blocks and slides

### Storage solutions



# Archives systems for blocks and slides

## Storage solutions




### Ergonomic design provides easy accessibility

Structure equipped with wheels to be conveniently moved in different environments.  
Or structure static equipped for transpallet use



### Easy to use, easy to work

Colored metal structure with white color. Drawers with color of your choice, and also available with closing doors with security lock (optional).



### High capacity

Detail of the archiving system for blocks with primary drawers and secondary drawers to maximize space.



### Safety first

Detail of the archiving system with drawers that can be completely removed.



### Half-capacity system

Storage for slides and blocks for half capacity, for lower volumes or to reach up to 2.5 m in height.  
But same characteristics.



### Security for your samples

Optional on request: archiving system equipped with safety and closing doors



### Technical support and after sale

Histo-Line offers a complete technical support for you. We find together the best solution for you and your instrument in certain and quick times.

Wide range of accessories and consumables on request.

### Histo-Line complete range of solutions for Archiving



Filing cabinet  
for blocks and slides



Filing cabinet  
for macro