



Innovative life science solutions

## Pre filled container Standard and Safe System



Disposable containers pre-filled with formaldehyde 4% based fixative, used in histology and cytology as primary fixative. The stability and the high quality of Histo-Line Laboratories formalin guarantee an highly effective and reproducible fixation.

**Wide range of available** sizes according to specimens to be fixed.

That containers are useful for the storage of pathology, industrial, and environmental samples, and these containers are made of **high-quality plastic** to resist breakage. Container lids screw on or snap into place easily and provide reliable **security for specimens** in transit or in storage.

### Formalin prefilled container safe system

Above the layer of formalin is an **inert liquid** and composed by non-hazardous reagents, maintains formalin features, protecting operators from formalin harmful vapors.

This product provides:

- fixed in formalin **without dangerous** leaks of formalin fumes
- **no changes** in lab
- **no changes in management** of specimens and protocols
- ready-to-use pre-filled containers
- **guarantee of safety** for operators

## Pre filled container Standard and Safe System

### Technical features

Code	Description	Packaging
07-124-01-S	Formalin prefilled container safe system (pressure cap)	22 pz. x 1000 ml
07-124-02-S	Formalin prefilled container safe system (screw cap)	224 pz x 20 ml.
07-124-03-S	Formalin prefilled container safe system (screw cap)	150 pz. x 30 ml
07-124-04-S	Formalin prefilled container safe system (screw cap)	150 pz. x 40 ml
07-124-06-S	Formalin prefilled container safe system (screw cap)	90 pz. x 60 ml
07-124-09-S	Formalin prefilled container safe system (screw cap)	90 pz. x 90 ml
07-124-12-S	Formalin prefilled container safe system (screw cap)	50 pz. x 120 ml
07-124-15-S	Formalin prefilled container safe system (screw cap)	50 pz. x 150 ml
07-124-16-S	Formalin prefilled container safe system (screw cap)	50 pz. x 160 ml
07-130-S	Containers transparent cap pre-filled safe system (pressure cap)	24 pz. x365 ml
07-131-S	Containers transparent cap pre-filled safe system (pressure cap)	22 pz. x1000 ml
07-133-S	Containers transparent cap pre-filled safe system (pressure cap)	16 pz. x500 ml
07-143-S	Containers transparent cap pre-filled safe system (pressure cap)	6 pz. x2000 ml
07-144-S	Containers transparent cap pre-filled safe system (pressure cap)	5 pz. x3000 ml
07-145-S	Containers transparent cap pre-filled safe system (pressure cap)	4 pz. x5000 ml
07-147-S	Containers transparent cap pre-filled safe system (pressure cap)	1 pz. x11.500 ml
07-148-S	Containers transparent cap pre-filled safe system (pressure cap)	1 pz. x26.500 ml

Code	Description	Packaging
07-124-12	Pre-filled containers with 10% neutral buffered formalin (screw cap) - 120 ml	50 pz. x 120 ml.
07-124-15	Pre-filled containers with 10% neutral buffered formalin (screw cap) - 150 ml	50 pz. x 150 ml.
07-124-16	Pre-filled containers with 10% neutral buffered formalin (screw cap) - 160 ml	50 pz. x 160 ml.
07-124-2	Pre-filled containers with 10% neutral buffered formalin (screw cap) - 20 ml	224 pz. x 20 ml.
07-124-3	Pre-filled containers with 10% neutral buffered formalin (screw cap) - 30 ml	150 pz. x 30 ml.
07-124-4	Pre-filled containers with 10% neutral buffered formalin (screw cap) - 40 ml	150 pz. x 40 ml.
07-124-6	Pre-filled containers with 10% neutral buffered formalin (screw cap) - 60 ml	90 pz. x 60 ml.
07-124-9	Pre-filled containers with 10% neutral buffered formalin (screw cap) - 90 ml	90 pz. x 90 ml.
07-130	Pre-filled containers with 10% neutral buffered formalin (pressure cap) - 365 ml	24 pz. x 365 ml.
07-131	Pre-filled containers with 10% neutral buffered formalin (pressure cap) - 1000 ml	22 pz. x 1000 ml.
07-133	Pre-filled containers with 10% neutral buffered formalin (pressure cap) - 500 ml	16 pz. x 500 ml.
07-143	Pre-filled containers with 10% neutral buffered formalin (pressure cap) - 2000 ml	6 pz. x 2000 ml.
07-144	Pre-filled containers with 10% neutral buffered formalin (pressure cap) - 3000 ml	5 pz. x 3000 ml.
07-145	Pre-filled containers with 10% neutral buffered formalin (pressure cap) - 5000 ml	4 pz. x 5000 ml.
07-147	Pre-filled containers with 10% neutral buffered formalin (pressure cap) - 11000 ml	1 pz. x 11500 ml.
07-148	Pre-filled containers with 10% neutral buffered formalin (pressure cap) - 26500 ml	1 pz. x26.500 ml



Viale G. di Vittorio, 30  
20090 Pantigliate (MI)  
ITALY  
Tel. +39 0255230061  
Fax +39 0255213764  
histoline@histoline.com  
[www.histoline.com](http://www.histoline.com)

Distributed by:

