

**Technical Data Sheet
Lab-O-Wax 52-54°C Paraffin for Histology**

CE IVD For In Vitro Diagnostic Use Manufacturer: Histo-Line Laboratories Srl

Code Product	Packaging
R0042 Paraffin Lab-O-Wax 52-54°C	1 Box with 6 bags of 2 kg
Code CND: W01030799	



Intended Use

Lab-O-Wax is a formulation of paraffin and plastic polymers intended for tissues embedding
 When melted, Lab-O-Wax produces a clear and bright solution, particle free
 Lab-O-Wax allows to cut very thin sections (1-2 micron) and excellent ribbons of sections.
 Lab-O-wax doesn't contain DMSO.
 For In Vitro Diagnostic Use

Physical Data

Fusion Point 52-54 °C
 Temperature Limit 60°C
 Appearance pellet

Validity

3 years

Storage Condition

Store at max. 45 °C. Store at max. 45 °C. At temperature higher than 45°C paraffin drops begin to aggregate.
 Keep the bags closed.

N.B. Although not reported in this technical sheet: compositions, precautionary measures, risk phrases, refer to storage, fire precautions, and more, refer to the Safety Data Sheet can be downloaded from our website: www.histoline.com or request it to: info@histoline.com

Milan Jul 2016