

**Technical Data Sheet  
Mono cassette**

**CE IVD** | *For In Vitro Diagnostic Use* | *Manufacturer: Histo-Line Laboratories Srl*

<b>Code Product</b>	<b>Packaging</b>
R4170 Mono Cassette White	1000 pcs.
R4157 Mono Cassette Blue	1000 pcs
R4154 Mono Cassette Yellow	1000 pcs
R4156 Mono Cassette Green	1000 pcs
R4155 Mono Cassette Pink	1000 pcs
R4159 Mono Cassette Red	1000 pcs
Codice CND: W05030501	



**Intended Use**

Cassettes for histological inclusions of biopsies and paraffin embedded samples, equipped with a double closing system, which guarantees a perfect sealing during processing.

**General features**

The casset presents a rough surface in the front portion to facilitate writing of graphite pencils. You can also use specific writing markers anti aggressive reagents and acids, or the use of special cassettes for imprint.

Internal dimensions: 30 x 25 x 5 mm, external dimensions: 40 x 28 x 6 mm.

These cassettes are compatibles with Histo-Line filing cabinet system and archives storage for high capacity.

**Functional features**

Material: acetal copolymer (POM polyacetal).

The material has excellent high temperature stability, to treatments in the microwave, solvents, acids and bases; then the material is optimal for the processing and the inclusion of pieces in paraffin.

**Validity**

Validity is not applicable to this type of product.

**Storage Condition**

The biocassettes not adopt special storage mode.

N.B. Although not reported in this technical sheet: compositions, precautionary measures, risk phrases, refer to storage, fire precautions, and more, refer to the Safety Data Sheet can be downloaded from our website: [www.histoline.com](http://www.histoline.com) or request it to: [info@histoline.com](mailto:info@histoline.com)

Milan Jul 2016